PhysTherm

Measure Temperature

Qosain Scientific's
PhysTherm measures
temperatures in a range as wide
as -260 °C to 750 °C with fast
response and notable accuracy.

It comes with a screw terminal connector that ensures support for a variety of K-type thermocouples. Simply hold the thermocouple probe in the environment of your choice, and using a guided software wizard, start visualizing and recording temperature in no time. PhysTherm utilizes an in-built cold junction compensation scheme which adds to the utility, making temperature measurement fast, reliable and as easy as the word "go"!



PhysLogger (needs to be purchased separately)

Typical Applications

Temperature measurements in any science lab or the outdoors

Monitoring temperature of perishable items

Temperature profling of the human body

Warming up of cryogenics

Modeling the behaviour of thermistor

The baffling Mpemba effect

Features

